



Crystallizing dishes


Code CNL001 - CNL009

Code CNN001 - CNN005

■ Endo Glassware provides with two models of cylindrical crystallizing dishes made of borosilicate or normal pressed glass that are ideal for drying or crystallizing wet solids thanks to their wide and plane base that allows a large evaporation surface.

■ Both models are available in several sizes to cover all users' necessities

DIMENSIONS - Crystallizing with spout



Borosilicate glass crystallizing dishes

Code CNL001 - CNL009



■ They can be heated and present a spout for an easy pouring.

Code	Capacity (approx.)	Diameter	Height
CNL001	40 mL	50 mm	30 mm
CNL002	60 mL	60 mm	35 mm
CNL003	100 mL	70 mm	40 mm
CNL004	150 mL	80 mm	45 mm
CNL005	250 mL	100 mm	50 mm
CNL006	500 mL	120 mm	65 mm
CNL007	1000 mL	150 mm	75 mm
CNL008	1500 mL	180 mm	90 mm
CNL009	3500 mL	230 mm	100 mm

Cristalizadores en vidrio prensado

Code CNN001 - CNN005



■ With reinforced rim to avoid chipping.

Code	Capacity (approx.)	Diameter	Height
CNN001	70 mL	80 mm	40 mm
CNN002	150 mL	100 mm	47 mm
CNN003	350 mL	120 mm	60 mm
CNN004	1000 mL	150 mm	70 mm
CNN005	1400 mL	190 mm	90 mm