

VITsafe™ safety wash bottles, narrow-mouth



Bottle made of PE-LD (or PP for acetone and MEK), PP spray tube.

More safety due to the durable safety imprint in accordance to Directive (EC) No.

1272/2008 (GHS), as well as with all important information:

- Material name in German, English, French and Spanish
- Chemical formula, CAS No., hazard pictogram, signal word
- Risk phrases (H phrases), safety phrases (P phrases), as well as NFPA Code

Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube. Practically no leakage or dripping with the VENT-CAP screw cap, the design of which prevents almost all static overpressure.



Imprint	Colour VENT-CAP	Volume ml	Thread Height*		Ø mm	PU	Cat. No.
			GL	mm			
Acetic acid	Red	500	25	180	74	6	1332979 / CKD537
Acetone	Red	250	25	135	65	12	1431829 / CKD543
Acetone	Red	500	25	180	74	12	1432829 / CKD544
Acetone	Red	1000	32	215	92	12	1433829 / CKB528
Acetonitrile	Red	500	25	180	74	6	1332969 / CKD536
Distilled Water	Blue	250	25	135	65	12	1331819 / CKD519
Distilled Water	Blue	500	25	180	74	12	1332819 / CKD524
Distilled Water	Red	1000	32	221	92	12	1333819 / CKD538
Ethanol	Orange	250	25	135	65	12	1331869 / CKD523
Ethanol	Orange	500	25	180	74	12	1332869 / CKD528
Ethanol	Red	1000	32	221	92	12	1333869 / CKD542
Ethyl acetate	Red	250	25	135	65	12	1331859 / CKD522
Ethyl acetate	Red	500	25	180	74	12	1332859 / CKD527
Ethyl acetate	Red	1000	32	221	92	12	1333859 / CKD541
Heptane	Red	500	25	180	74	6	1332899 / CKD531
i-Hexane	Red	500	25	180	74	6	1332909 / CKD532
Isopropanol	Yellow	250	25	135	65	12	1331849 / CKD521
Isopropanol	Yellow	500	25	180	74	12	1332849 / CKD526
Isopropanol	Red	1000	32	221	92	12	1333849 / CKD540
Methanol	Green	250	25	135	65	12	1331839 / CKD520
Methanol	Green	500	25	180	74	12	1332839 / CKD525
Methanol	Red	1000	32	221	92	12	1333839 / CKD539
Methyl ethyl ketone (MEK)	Red	500	25	180	74	6	1432989 / CKD545
Methylene chloride	Transparent	500	25	180	74	6	1332879 / CKD529
N,N-Dimethylformamide	Red	500	25	180	74	6	1332889 / CKD530
Pentane	Red	500	25	180	74	6	1433959
Tetrahydrofuran (THF)	Red	500	25	180	74	6	1332939 / CKD533
Toluene	Red	500	25	180	74	6	1332949 / CKD534
Xylene	Red	500	25	180	74	6	1332959 / CKD535

* Height without spray tube

Other variations available upon request.



VITsafe™ safety wash bottles, wide-mouth



Bottle made of PE-LD (or PP for acetone and MEK), PP spray tube.

More safety due to the durable safety imprint in accordance to Directive (EC) No. 1272/2008 (GHS), as well as with all important information:

- Material name in German, English, French and Spanish
- chemical formula, CAS No., hazard pictogram, signal word
- Risk phrases (H phrases), safety phrases (P phrases), as well as NFPA Code

Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube. Practically no leakage or dripping with the VENT-CAP screw cap, the design of which prevents almost all static overpressure.

Imprint	Colour VENT-CAP	Volume ml	Thread GL	Height* mm	Ø mm	PU	Cat. No.
Acetic acid	Red	500	45	166	76	6	1352979 / CKB522
Acetone	Red	250	45	146	58	12	1451829/CKB529
Acetone	Red	500	45	166	76	12	1452829/CKB530
Acetone	Red	1000	63	226	91	12	1453829/CKB532
Acetonitrile	Red	500	45	166	76	6	1352969/CKB521
Distilled Water	Blue	250	45	146	58	12	1351819/CKB504
Distilled Water	Blue	500	45	166	76	12	1352819/CKB509
Distilled Water	Red	1000	63	226	91	12	1353819/CKB523
Ethanol	Orange	250	45	146	58	12	1351869/CKB508
Ethanol	Orange	500	45	166	76	12	1352869/CKB513
Ethanol	Red	1000	63	226	91	12	1353869/CKB527
Ethyl acetate	Red	250	45	146	58	12	1351859/CKB507
Ethyl acetate	Red	500	45	166	76	12	1352859/CKB512
Ethyl acetate	Red	1000	63	226	91	12	1353859/CKB526
Heptane	Red	500	45	166	76	6	1352899/CKB516
i-Hexane	Red	500	45	166	76	6	1352909/CKB517
Isopropanol	Yellow	250	45	146	58	12	1351849/CKB506
Isopropanol	Yellow	500	45	166	76	12	1352849/CKB511
Isopropanol	Red	1000	63	226	91	12	1353849/CKB525
Methanol	Green	250	45	146	58	12	1351839/CKB505
Methanol	Green	500	45	166	76	12	1352839/CKB510
Methanol	Red	1000	63	226	91	12	1353839/CKB524
Methyl ethyl ketone (MEK)	Red	500	45	166	76	6	1452989/CKB531
Methylene chloride	Transparent	500	45	166	76	6	1352879/CKB514
N,N-Dimethylformamide	Red	500	45	166	76	6	1352889/CKB515
Pentane	Red	500	45	166	76	6	1453959/CKD546
Tetrahydrofuran (THF)	Red	500	45	166	76	6	1352939/CKB518
Toluene	Red	500	45	166	76	6	1352949/CKB519
Xylene	Red	500	45	166	76	6	1352959/CKB520

* Height without spray tube

Other variations available upon request.


VENT-CAP wash bottle caps, PP

Screw cap and spray tube from PP.

Practically no leakage or dripping with the VENT-CAP screw cap, the design of which prevents almost all static overpressure.

Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube.

For bottle size ml	Thread GL	Colour	PU	Cat. No.
250	25	Red	12	834111
250	25	Blue	12	834121
250	25	Green	12	834131
250	25	Orange	12	834141
250	25	Yellow	12	834151
250	45	Red	12	834311
250	45	Blue	12	834321
250	45	Green	12	834331
250	45	Orange	12	834341
250	45	Yellow	12	834351
500	25	Transparent	12	834102
500	25	Red	12	834112
500	25	Blue	12	834122
500	25	Green	12	834132
500	25	Orange	12	834142
500	25	Yellow	12	834152
500	45	Transparent	12	834302
500	45	Red	12	834312
500	45	Blue	12	834322
500	45	Green	12	834332
500	45	Orange	12	834342
500	45	Yellow	12	834352
1000	32	Red	12	834213
1000	63	Red	12	834413





Wash bottles with imprint, PE-LD/PP



Narrow- / wide-mouth bottles made of PE-LD, transparent.

Screw cap and spray tube from PP.

Imprinted with "Distilled Water" in German, English, French and Spanish.

Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube.

Upon request, also with other imprints for non-hazardous media according to the REACH Directive.

Volume ml	Thread GL	Height* mm	Ø mm	PU	Cat. No.
250	25	135	65	12	133181/CKD512
250	45	146	58	12	135181/CKB504
500	25	180	74	12	133281/CKD513
500	45	166	76	12	135281/CKB502
1000	32	221	92	12	133381/CKD514
1000	63	226	91	12	135381/CKB503

* Height without spray tube



Wash-bottles, PP



Narrow- / wide-mouth bottles made of PP, transparent. Screw cap and spray tube from PP.

Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube.

Volume ml	Thread GL	Height* mm	Ø mm	PU	Cat. No.
250	25	135	65	12	94993/CKD550
250	45	146	58	12	93793/CKB534
500	25	180	74	12	95093/CKD552
500	45	166	76	12	93993/CKB536
1000	32	215	92	12	95193/CKD554
1000	63	226	91	12	94193/CKB538

* Height without spray tube

Wash-bottles, PE-LD/PP



Narrow- / wide-mouth bottles made of PE-LD, transparent. Screw cap and spray tube from PP.

Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube.

Volume ml	Thread GL	Height* mm	Ø mm	PU	Cat. No.
50	18	85	37	24	94588/CKD547
100	18	114	43	24	94688/CKD548
250	25	135	65	12	94988/CKD549
250	45	146	58	12	93788/CKB533
500	25	180	74	12	95088/CKD551
500	45	166	76	12	93988/CKB535
1000	32	221	92	12	95188/CKD553
1000	63	226	91	12	94188/CKB537

* Height without spray tube



Wash-bottles, coloured, PE-LD/PP

Narrow-mouth bottles, made from PE-LD. Available in four different colours to facilitate ready identification. Screw cap and spray tube from PP.

Precise spray jet and optimised medium backflow through the continuous and finely drawn tip of the spray tube.

Colour	Volume ml	Thread GL	Height* mm	Ø mm	PU	Cat. No.
red	250	25	135	65	5	132603/CKD504
red	500	25	180	74	5	132703/CKD508
red	1000	32	221	92	5	132803/CKD560
green	250	25	135	65	5	132605/CKD505
green	500	25	180	74	5	132705/CKD509
green	1000	32	221	92	5	132805/CKD561
yellow	250	25	135	65	5	132606/CKD506
yellow	500	25	180	74	5	132706/CKD510
yellow	1000	32	221	92	5	132806
blue	250	25	135	65	5	132608/CKD507
blue	500	25	180	74	5	132708/CKD511
blue	1000	32	221	92	5	132808/CKD563
Set: red, green, yellow, blue (1 item each)	500	25	180	74	1	1327111/CKD555
Set: red, green, yellow, blue (1 item each)	1000	32	221	92	1	1328111

* Height without spray tube






Wash-bottles for individual labeling, with VENT-CAP, PE-LD/PP

Narrow- and wide-mouth bottles made from PE-LD, transparent. VENT-CAP screw cap and spray tube made from PP.

Practically no leakage or dripping of the liquid due to the VENT-CAP screw cap, that prevents almost all build-up of static overpressure. The spray tube enables a precise spray jet and optimised medium backflow.

Volume ml	Thread GL	Height* mm	Ø mm	PU	Cat.No.
500	25	180	74	6	133211
500	45	166	76	6	135211

*Height without spray tube



Wash bottle caps, PP

Screw cap and spray tube with drawn-out tip, made of PP. With precise spray jet and optimised medium backflow.

Please note the different lengths of the spray tubes (short - medium - long):

- Cat. No. 833001 – short
- Cat. No. 833101 – short
- Cat. No. 833102 – medium
- Cat. No. 833203 – long
- Cat. No. 833301 – short
- Cat. No. 833302 – medium
- Cat. No. 833403 – long

Thread GL	PU	Cat.No.
18	24	833001
25	12	833101
25	12	833102
32	12	833203
45	12	833301
45	12	833302
63	12	833403