

Reagent bottles, PP



Transparent.

With NS stopper made of PP.

Stopper type A: With square-knob cap and red core.

Stopper type B: With octagonal-knob cap and red core.

Good chemical resistance, ideal for long-term storage of liquids.

Autoclavable at 121 °C (2 bar) according to DIN EN 285.



Volume ml	Neck NS	Height mm	Ø mm	Stopper	PU	Cat. No.
100	14/23	106	52	A	20	100394/CHJ501
100	29/32	111	55	B	20	101594/CHD501
250	19/26	138	70	A	20	100494
250	34/35	144	73	B	20	101694
500	24/29	172	87	A	10	100594/CHJ503
500	45/40	183	87	B	10	101794/CHD503
1000	29/32	213	108	A	10	100694/CHJ504
1000	60/46	214	108	B	10	101894/CHD504

Standard joint stoppers, PP



Stopper type A: With square-knob cap and red core.

Stopper type B: With octagonal-knob cap and red core.

Autoclavable at 121 °C (2 bar) according to DIN EN 285.



NS	Model	PU	Cat. No.
10/19	A	10	90694/BMB500
12/21	A	10	90794/BMB501
14/23	A	10	90894/BMB502
19/26	A	10	90994/BMB503
24/29	A	10	91094/BMB504
29/32	A	10	91194/BMB505
29/32	B	10	92194/BMB509
34/35	B	10	91294/BMB506
45/40	B	10	91394/BMB507
60/46	B	10	91494/BMB508