

## BURETTE FLASK

clear


## BÜRETENFLASCHE

klar

**Cat. No. 2100 B**

| Catalog No.<br>Katalog-Nr. | Capacity (ml)<br>Kapazität | d (mm) | h (mm) | SJ<br>NS | Packing (pcs)<br>Packung (Stück) |
|----------------------------|----------------------------|--------|--------|----------|----------------------------------|
| N555436801046              | 2000                       | 160    | 200    | 29/32    | 6                                |

CSN  
70 4118



## MEASURE GRADUATED, BELL SHAPED

with spout, graduated

## MENSUREN GRADUIERT, GLOCKENFÖRMIG

mit Ausguss

**Cat. No. 1686**

| Catalog No.<br>Katalog-Nr. | Capacity (ml)<br>Kapazität | d1 (mm) | d2 (mm) | h (mm) | Packing (pcs)<br>Packung (Stück) |
|----------------------------|----------------------------|---------|---------|--------|----------------------------------|
| 1632432540818              | 5                          | 26      | 22      | 70     | 10                               |
| 1632432540919              | 10                         | 30      | 25      | 75     | 10                               |
| 1632432541123              | 25                         | 39      | 30      | 90     | 10                               |
| 1632432541125              | 50                         | 50      | 38      | 110    | 10                               |
| 1632432541230              | 100                        | 62      | 49      | 130    | 24                               |
| 1632432541338              | 250                        | 84      | 67      | 160    | 10                               |
| 1632432541443              | 500                        | 102     | 80      | 200    | 10                               |
| 1632432541644              | 1000                       | 126     | 100     | 240    | 2                                |
| 1632432541746              | 2000                       | 148     | 116     | 310    | 8                                |

