

## Spare parts and bottles



### Standard and telescopic feed tubings

Supplied either cut to standard size or uncut, by the meter.

Optional telescopic tubing for instant fitting to bottle size.

Description	Length	int. Ø	Cat. No.
-------------	--------	--------	----------

#### Replacement feed tubing, PTFE

Fits Calibrex™ 520/525/530 up to 10 mL, cut	300 mm	5 mm	511.707
Fits Calibrex™ 520/525/530 up to 10 mL, by the meter, uncut	specify	5 mm	511.709

#### Replacement feed tubing, FEP

For Calibrex™ 525/530 as of 25 mL, cut	350 mm	7 mm	525.350
For Calibrex™ 525/530 as of 25 mL, by the meter, uncut	specify	7 mm	525.706

#### Telescopic feed tubing, FEP

Fits Calibrex™ 520/525/530 up to 10 mL	150 - 255 mm	1 / pk	1.525.352
Fits Calibrex™ 525/530 as of 25 mL	195 - 345 mm	1 / pk	1.525.355



### Alternative setting button

Classical screw-type button easily replaces the original sliding cursor, according to user preference.

### Calibration seal stickers

Autoclavable seal sticker protects access to the calibration mechanism of Calibrex™ 525/530 models.

Description	Colour	Packaging	Cat. No.
-------------	--------	-----------	----------

#### Screw type volume setting button

Fits Calibrex™ 525 dispensers	Yellow	1 / pk	1.525.918
Fits Calibrex™ 530 dispensers	Red	1 / pk	1.530.918

#### Safety seal sticker

Fits Calibrex™ 525/530 dispensers		25 / pk	1.525.525
Fits Calibrex™ 525/530 dispensers		100 / pk	1.525.526



### Bottle-neck adapters

Each Calibrex™ dispenser comes along with three bottle-neck adapters. To fit additional bottle-neck diameter, order the appropriate adapter size or combine two adapters to reach the adequate solution.

Type	Ext. Ø of bottle neck	Cat. No.
------	-----------------------	----------

#### Adapters for Calibrex™ organo 525 and solutae 530

Threaded, PP	32 - 25 mm	1.525.GL25
Threaded, PP	32 - 28 mm	1.525.GL28
Threaded, PP	45 - 32 mm	1.525.GL32
Threaded, PP	45 - 38 mm	1.525.GL38
Threaded, PP	45 - 40 mm	1.525.GL40

#### Adapters for Calibrex™ universal 520

Threaded, PP	Ø 22, 25, 28, 30, 34, 36, 38, 40, 45 mm	GLP + Ø
Threaded, PTFE	Ø 36, 38, 40, 45 mm	GLT + Ø



### Reagent bottles

Glass and polyethylene reagent bottles supplied with its PP screw cap. Suitable for all bottle top dispensers.

Corresponding neck adapters supplied with the dispensers.

Shape	Volume	Neck, ext. Ø	Cat. No.
-------	--------	--------------	----------

#### Amber glass<sup>1</sup>

Square	100 mL	32 mm	314.0100
Square	250 mL	32 mm	314.0250
Square	500 mL	32 mm	314.0500
Square	1000 mL	45 mm	314.1000
Round	2500 mL	45 mm	314.2500

#### Amber glass<sup>2</sup>, with handle

Round	2500 mL	45 mm	314.2500H
-------	---------	-------	-----------

#### PE coated amber glass<sup>2</sup>

Square	500 mL	32 mm	314.0500PE
Square	1000 mL	45 mm	314.1000PE
Round	2500 mL	45 mm	314.2500PE

Shape	Volume	Neck, ext. Ø	Cat. No.
-------	--------	--------------	----------

#### Clear borosilicate glass<sup>1</sup>

Round	250 mL	45 mm	314.0250C
Round	500 mL	45 mm	314.0500C
Round	1000 mL	45 mm	314.1000C
Round	2000 mL	45 mm	314.2000C

#### Polyethylene<sup>2</sup>

Square	250 mL	25 mm	315.0250
Square	500 mL	25 mm	315.0500
Square	1000 mL	32 mm	315.1000
Square	2500 mL	45 mm	315.2500

<sup>1</sup> autoclavable <sup>2</sup> not autoclavable