



Digital vortex multi-tube mixer

Code LGB004

The Nahita Blue digital vortex multi-tube mixer adopts a simple programming design with the advantage of being able to mix and process up to 54 samples at a time, saving time and effort.

- Maintenance-free brushless motor.
- Equipped with emergency stop button, to immediately stop operation of the equipment in case of need.
- Smooth start with uniform acceleration
- It incorporates 10 storage programmes.
- It incorporates the 'Pulse' function, which allows you to select the operating and stop time. It also has the option of setting the number of cycle repetitions.
- 3-inch LCD display, with online settings and clear operating parameters.
- It incorporates an acoustic warning at the end of the cycle.
- Graduated rods with a scale from 0 to 15 cm for correct alignment.
- Made of metal, with four suction cup rubber feet.



Technical specifications

Code	LGB004
Speed range	200 - 2500 rpm (10rpm increments)
Speed accuracy	±15 rpm
Time range	1-59 s / 1-999 min
Oscillation mode	Circumference
Circumferential diameter	3.6 mm
Display	Backlit LCD
Programmes	0-9
Maximum load	5 kg
Input power supply	100/240 V; 50/60 Hz
Power	90 W
Level of protection	IP20
Weight	24 Kg
Dimensions (W x D x H)	425x300x445 mm

Built-in accessories



Description	1 sponge rack for 54 tubes 1.5/2 mL	1 sponge rack for 54 x 15 mL tubes	1 sponge rack for 20 x 50 mL tubes
--------------------	-------------------------------------	------------------------------------	------------------------------------