



## Draining rack

Code ZCZ001, ZCZ002, ZCZ003, ZCZ016

As essential elements in the daily lab work, table and wall draining racks are used to drain and dry all kind of laboratory material (glassware, plasticware, etc.) of small and medium size (beakers, measuring cylinders, flasks, tubes, etc.). Available both table and wall models.

### Table models

Code ZCZ001, ZCZ002, ZCZ016



- Made of plastified wire, they are supplied with silicone rubber protectors and extractable tray in the base.
- They can stand on a table surface or can be hanged on the wall.

Code	Positions	Dimensions	Type
ZCZ001	32	470x360x160 mm	Table/Wall
ZCZ002	55	640x360x140 mm	Table/Wall

#### Accessories

Code	Description
ZCZ016	Tray for draining rack ZCZ001, ZCZ002

### Wall model

Code ZCZ003



- Made of polypropylene, pins can be easily detached allowing adapting to material with special size or shape.
- It is provided with drainpipe for emptying the lower tray.

Code	Positions	Dimensions	Type
ZCZ003	27	550x400x120 mm	Wall