




Spherical flasks

Code CDG001_007, CDJ001_007

- We offer this wide range of spherical flasks made of thermal shock resistant borosilicate glass.
- With narrow neck and a white marking area, these flasks are ideal for distillation processes or those reactions in which a uniform sample heating is necessary
- They are available in different capacities to cover all users' necessities.

FLASK, glass

d1: diameter
d2: diameter
h: height



Flat bottom

Code CDG001_007



Technical specifications

Code	Cap. mL	d1. mm	d2. mm	h. mm
CDG001	50	50	22	105
CDG002	100	65	22	108
CDG003	250	85	34	140
CDG004	500	105	34	180
CDG005	1000	135	43	210
CDG006	2000	160	52	250
CDG007	3000	188	59	282

Round bottom

Code CDJ001_007



Technical specifications

Code	Cap. mL	d1. mm	d2. mm	h. mm
CDJ001	50	51	22	105
CDJ002	100	64	22	115
CDJ003	250	85	34	145
CDJ004	500	105	34	175
CDJ005	1000	131	42	210
CDJ006	2000	166	50	260
CDJ007	5000	223	50	345