



Designed to hold, serve and carry instruments and consumables, it is perfectly adapted to be used in ward, nursing, medication, or any department in the hospital matching infection control requirements.

Sturdy and easy maneuverable it has a loading capacity of 60 kg. Easy for cleaning and disinfection resistant to major cleaning solutions used in hospital. This trolley is designed and produced in France thanks to our 60 years expertise knowledge.



DURABILITY

QUALITY

STURDINESS

- Easy to be exported worldwide due to the design and flat packing for an optimal freight cost.
- Customisable and modular with wide range of optionnal accessories that are attachable even after installation in hospital.
- Quick assembly by technical staff in the hospital
- Guarantee : 2 years



AGENCINOX

6 rue Marcel Brault - ZI
89110 Aillant sur Tholon
FRANCE

Tel. +33 (0)3 86 92 39 27
Fax. + 33 (0)3 86 92 39 07

Mail : contact@agencinox.com

site : www.agencinox.com



DRESSING TROLLEYS

Technical specifications :

- Stamped shelves 750 x 500 mm S.S 18/10 AISI 304
- S.S 18/10 AISI 304 Tube Ø 25 mm
- Shelves height from floor : 300 and 820 mm for 2 shelves
- Shelves height from floor : 200, 510 and 820 mm for 3 shelves
- Wheels Ø 75 mm
- 4 Bumpers 95 mm
- Assembly with S.S screws and aluminium rods
- Delivered dismantled with assembling instructions and tools
- CE medical device class I
- Packing : 62 x 11 x 103 cm

Accessories :



Swivel plastic bag holder 10L
Ref : 114



Sharps container holder

Ref : 123 (without box)

**Small drawer
255 x 400 x 115 mm**
Ref : 153



**Large drawer
550 x 400 x 115 mm**
Ref : 154

Guard rail (length)

Ref : 741E6

Guard rail (width)

Ref : 741E5

Clip for guard rail

Ref : P27

REFERENCES	DESCRIPTIONS	DIMENSIONS	WEIGHT
741	Dressing trolley 2 shelves	750 x 500 mm	11.5 kg
745	Dressing trolley 2 shelves with guard rails	750 x 500 mm	12.5 kg
746	Dressing trolley 3 shelves	750 x 500 mm	15 kg
747	Dressing trolley 3 shelves with guard rails	750 x 500 mm	17.5 kg

Options :

- Ref : X ■ Ref : A ■ Ref : F
- Wheels Ø 100mm Antistatic wheels 2 wheels with brakes

Example : trolley 750 x 500 mm with antistatic wheels + brakes = Ref : **741AF**



AGENCINOX

6 rue Marcel Brault - ZI
89110 Aillant sur Tholon
FRANCE

Tel. +33 (0)3 86 92 39 27
Fax. + 33 (0)3 86 92 39 07

Mail : contact@agencinox.com

site : www.agencinox.com