



## Quartz crucible without lid

Code CPH016\_ CPH019



- Made of quartz gel and powder.
- Heat tolerance up to 1200 °C.

Code	Capacity mL	ø diam.x height	Type
CPH016	15	38x28	Without lid
CPH017	25	44x36	Without lid
CPH018	50	53x45	Without lid
CPH019	100	67x55	Without lid

## Quartz crucible with lid

Code CPH020\_ CPH023



- Made of quartz gel and powder.
- Heat tolerance up to 1200 °C.

Code	Capacity mL	ø diam.x height	Type
CPH020	15	38x26	With lid
CPH021	25	44x36	With lid
CPH022	50	53x45	With lid
CPH023	100	67x55	With lid