



Narrow-mouth bottles, PP



Transparent. With high shoulders.
 With screw cap made of PP.
 Good chemical resistance, ideal for long-term storage of liquids.
 Autoclavable at 121 °C (2 bar) according to DIN EN 285.
 Suitable for contact with foodstuffs according to regulation (EU) No. 10/2011.

| Volume ml | Thread GL | Height mm | Ø mm | PU | Cat. No. |
|-----------|-----------|-----------|------|----|----------|
| 250 | 25 | 135 | 65 | 12 | 94994 |
| 500 | 25 | 180 | 74 | 12 | 95094 |
| 1000 | 32 | 215 | 92 | 12 | 95194 |



Wide-mouth bottles, PP



Transparent.
 With screw cap made of PP.
 Good chemical resistance, ideal for long-term storage of liquids.
 Autoclavable at 121 °C (2 bar) according to DIN EN 285.
 Wider mouth for simple filling, also ideal for powders and paste-like materials.
 Suitable for contact with foodstuffs according to regulation (EU) No. 10/2011.

| Volume ml | Thread GL | Height mm | Ø mm | PU | Cat. No. |
|-----------|-----------|-----------|------|----|----------|
| 250 | 45 | 146 | 58 | 12 | 93794 |
| 500 | 45 | 166 | 76 | 12 | 93994 |
| 1000 | 63 | 226 | 91 | 12 | 94194 |