




PP tube racks - for tubes


Code **BJB078, BJB079, BJB080, BJB081, BJB082, BJB083, BJB084, BJB085, BJB086, BJB087, BJB088, BJB089, BJB090, BJB091, BJB092, BJB093, BJB094, BJB095, BJB096**



- Made of high-quality polypropylene
- Resistant to alcohols and mild organic solvents
- Stable temperature ranges from -196 °C to 121 °C
- Suitable for use in water baths as it remains submerged
- Supplied flat for easy assembly

Technical specifications



Reference	Tubes	Alignment	Dimensions (mm)
BJB078	50 mL x 20 / 15 mL x 30	5 x 4 / 6 x 5	208 x 172 x 60



- Made of high-quality polypropylene
- Resistant to alcohols and mild organic solvents
- Stable temperature ranges from -196 °C to 121 °C
- Alphanumeric markings on the holes
- Suitable for use in water baths as it remains submerged
- Supplied flat for easy assembly


Technical specifications

Reference	Tubes	Alignment	Dimensions (mm)
BJB081	∅ 13 mm x 90	15 x 6	247 x 107 x 60
BJB082	∅ 13 mm x 90	15 x 6	247 x 107 x 60
BJB083	∅ 13 mm x 90	15 x 6	247 x 107 x 60
BJB084	∅ 13 mm x 90	15 x 6	247 x 107 x 60
BJB085	∅ 16 mm x 60	12 x 5	247 x 107 x 60
BJB086	∅ 16 mm x 60	12 x 5	247 x 107 x 60
BJB087	∅ 16 mm x 60	12 x 5	247 x 107 x 60
BJB088	∅ 16 mm x 60	12 x 5	247 x 107 x 60
BJB089	∅ 21 mm x 40	10 x 4	247 x 107 x 60
BJB090	∅ 21 mm x 40	10 x 4	247 x 107 x 60
BJB091	∅ 21 mm x 40	10 x 4	247 x 107 x 60
BJB092	∅ 21 mm x 40	10 x 4	247 x 107 x 60
BJB093	∅ 30 mm x 21	7 x 3	247 x 107 x 60
BJB094	∅ 30 mm x 21	7 x 3	247 x 107 x 60
BJB095	∅ 30 mm x 21	7 x 3	247 x 107 x 60
BJB096	∅ 30 mm x 21	7 x 3	247 x 107 x 60





PP tube racks - for tubes




- Made of high-quality polypropylene
- Resistant to alcohols and mild organic solvents
- Stable temperature ranges from -196 °C to 121 °C
- Suitable for use in water baths as it remains submerged
- Supplied flat for easy assembly

Technical specifications

Reference	Tubes	Alignment	Dimensions (mm)	
BJB079	15 mL x 50	10 x 5	207 x 130 x 56	
BJB080	50 mL x 25	5 x 5	201 x 171 x 60	

PP tube racks - for microtubes

Code [BJL018](#), [BJL019](#), [BJL020](#), [BJL021](#), [BJL022](#), [BJL023](#), [BJL024](#), [BJL025](#), [BJL026](#)



- Made of high-quality polypropylene
- Resistant to alcohols and mild organic solvents
- Stable temperature ranges from -196 °C to 121 °C
- Removable lid includes a writable area for inventory lists

Technical specifications

Reference	Tubes	Alignment	Dimensions (mm)	
BJL018	0,2 mL x 96	12 x 8	120 x 82 x 30	
BJL019	0,5 / 1,5 mL x 60	12 x 5	198 x 73 x 49	
BJL020	0,5 / 1,5 mL x 60	12 x 5	198 x 73 x 49	
BJL021	0,5 / 1,5 mL x 60	12 x 5	198 x 73 x 49	
BJL022	0,5 / 1,5 mL x 60	12 x 5	198 x 73 x 49	
BJL023	0,5 / 1,5 mL x 96	12 x 8	198 x 112 x 49	
BJL024	0,5 / 1,5 mL x 96	12 x 8	198 x 112 x 49	
BJL025	0,5 / 1,5 mL x 96	12 x 8	198 x 112 x 49	
BJL026	0,5 / 1,5 mL x 96	12 x 8	198 x 112 x 49	