



CULTURE MULTIWELL PLATES, TREATED, STERILE

Code [BRL001](#) [BRL002](#), [BRL003](#), [BRL004](#), [BRL005](#), [BRL006](#), [BRL007](#)

Sterile Endo Plasticware multiwell culture plates are essential consumables in biology, microbiology, and biotechnology laboratories. They are specially designed for the simultaneous cultivation and analysis of multiple samples, enabling efficient, reproducible, and easily comparable assays in a compact and versatile format.

- Made of 100% virgin polystyrene, transparent colour
- Treated Surface
- Uniform volume and growth surface to maximize the usable growth area in all Wells
- Alphanumeric coding for easy identification
- Electron beam (E-beam) sterilized
- DNase, RNase, and pyrogen free
- Lids included
- Individually packaged



Technical Specifications

Code	Wells	Bottom	Working Volume (ml)	Growth Area (cm ²)
BRL001	6	Flat	1.90-2.90	9.60
BRL002	12	Flat	0.76-1.14	3.85
BRL003	24	Flat	0.38-0.57	1.93
BRL004	48	Flat	0.19-0.29	0.83
BRL005	96	Flat	0.075-0.20	0.33
BRL006	96	U-shaped	0.075-0.20	0.58
BRL007	96	V-shaped	0.075-0.20	