



JIZERSKÁ PORCELÁNKA, s.r.o.

Producer of laboratory and technical porcelainsince1947




PRODUCT SPECIFICATIONS SHEET

Casserole with porcelain handle

Jizerská porcelánka s.r.o. certifies, that the specified products meet all specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

Sign	JIPO Code	V (ml)	H (mm)	D (mm)
209/1	641 322 209 101	22	26	45
209/2	641 322 209 201	84	39	67
209/3	641 322 209 301	176	49	85
209/4	641 322 209 401	270	51	101
209/5	641 322 209 501	450	64	120
209/6	641 322 209 601	960	86	150



- Intended uses:** Laboratory applications. Can be used in microwave, dishwasher, direct fire only with ceramic wire flame diffuser. Do not use for induction heating.
- Material:** Porcelain C110
- Shape:** Flat bottom,withspoutand handle, without graduation



JIZERSKÁ PORCELÁNKA, s.r.o.

Producer of laboratory and technical porcelainsince1947



PRODUCT SPECIFICATIONS SHEET

Casserole with porcelain handle

Compliant to: ISO material for the manufacture lies within the group C110 according to ČSN EN 60672 –3
Product dimension tolerances are controlled by DIN 40680
Acid resistance acc. to ISO 1775

Properties: Unglazed products can be used at temperatures up to 1350°C / 2462°F and glazed ones at temperatures of up to 1000°C / 1,832°F (Δt 230°C)
Products have high resistance to acids and bases
Density 2.4 g/cm³
Bending strength of unglazed products (70MPa)
Thermal conductivity (W.m-1.K-1) = 1.3

The compliance of demands is proved by regular tests and followed noticed tests from an independent laboratory:

Test according ČSN EN60672-2
Determination of thermal shock resistance Δt 230°C
Test in liquid nitrogen
Test resistance to acids and bases

Disclaimer: *Jizerská porcelánka s.r.o. wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use. If there's any doubt ask us.*