

TAPERED STOPPERS - LONG SERIES (slight taper)

Versilic grade:

Made of pure silicone. Taste- and odor-free. Amber-colored.

Excellent chemical inertia, biologically neutral. Resistance to aging, ozone and repeated sterilization.

Temperature resistance from -60°C to 200°C (peaks acceptable up to 220°C).




Iso-Versinic grade:

Suitable for use in temperatures from -20°C to 200°C (peaks acceptable up to 300°C).

Exceptional chemical stability. Resistance to numerous acids and aggressive chemicals, even at high concentrations and temperatures, notably hydrocarbons, oils and strong acids, aromatic and chlorinated solvents.

Stopper geometry varies according to the item number. Hollow-bottom design makes for better seal:

Type	Dimensions (mm)			VERSILIC			ISO VERSINIC		
	Dia. base	Dia. top	Height	Pkg. Qty.*	Item no.	Shape	Pkg. Qty.*	Item no.	Shape
000 0000	3	7	19	500	261002	F1	50	271000	F1
000 000	3.7	6	5	500	261003	F1	50	271002	F1
000 00	5	7	8	500	261004	F1			
000 0	7	10	9	500	261006	F1	50	271006	F1
000	8.5	10	13	500	261008	F1	50	271008	F1
000 sharp	9.2	11.2	13.2						
00	6	10	17.5	500	261014	F1	50	271014	F1
0	9	12	19	500	261016	F1	50	271016	F1
1 slight	10	14	19	500	261020	F1	25	271020	F1
1	12	15	19	500	261024	F1	25	271024	F1
1 n	11	16	25	500	261028	F1			
2 slight	13	16	21	250	261030	F1	25	271030	F1
2	14	17	21	250	261034	F1	25	271034	F1
2 n	12	18	25	250	261040	F1			
3	16	19	23	100	262042	F2	10	272042	F2
4	18	21	25	100	262050	F2	10	272050	F2
5	19	23	28	50	262058	F2	10	272058	F2
6	22	26	30	50	262066	F2	10	272066	F2
7	24	29	32	50	262074	F2	10	272074	F3
8	26	32	34	25	262082	F2	5	272082	F3
9	28	34	36	25	262090	F3	5	272090	F3
10	31	37	41	25	262098	F3	5	272098	F3
11	34	40	43	25	262106	F3	5	272106	F3
12	36	43	45	25	262114	F3	5	272114	F3
13	40	49	49	20	262122	F3	5	272122	F3
14	46	55	51	10	262130	F3	5	272130	F3
15	50	61	56	10	262138	F3	5	272138	F3
16	54	63	61	10	262146	F3			
17	61	71	65	5	262154	F3			
18	69	84	65	5	262162	F3			
19	80	100	67	5	262170	F3			
20 slight	93	112	88						
20	100	130	85	1	262186	F3			
21	120	170	130	1	262194	F3			



These items are not in stock and must be specially ordered.

Pkg. Qty. = Delivered in packs of items.