



CHEMICAL SCOOP

Stainless steel 1.4301

Content in cm ³	120	450	1000	1500	2500
Ø in mm	50	75	100	120	150
Scoop length in mm	100	150	200	230	250
Handle length in mm	95	105	105	105	105
	80.758.050	80.758.075	80.758.100	80.758.120	80.758.150

CHEMICAL SCOOP

Aluminium

Length in mm	Content in ml	
160	75	80.748.160
185	113	80.748.185
210	210	80.748.210
250	300	80.748.250
265	410	80.748.265
310	650	80.748.310
350	1100	80.748.350
400	1870	80.748.400
430	2580	80.748.430



GRADUATED MEASURING JUG

Stainless steel 1.4301, heavy version

Content in ml	1000	1500	2000
Height in mm	140	170	170
Ø in mm	100	110	130
	80.718.100	80.718.150	80.718.200



BUCKET

Stainless steel 1.4301, graduated, with flat bottom

Content in l	6	10	12	15
bucket	80.730.006	80.730.010	80.730.012	80.730.015
cover	80.731.006	80.731.010	80.731.012	80.731.015

