


Magnetic stirring bar retrievers, flexible, PTFE



Flexible magnetic stir bar retriever with a total length of 330 mm. Magnet encapsulated. $\emptyset \times L$: 12.5 x 51 mm. Due to the high flexibility, the magnet stir bar can be retrieved from inaccessible locations, e.g., from the water trap in a laboratory sink. High chemical resistance, simple to clean.

| Length mm | PU | Cat. No. |
|--------------|----|----------|
| 330 | 1 | 318597 |



Magnetic stirring bar retrievers, PTFE



With PTFE encapsulated magnetic core. Straight shape. High chemical resistance, simple to clean.

| Length mm | PU | Cat. No. |
|--------------|----|----------|
| 150 | 1 | 122097 |
| 250 | 1 | 122197 |
| 350 | 1 | 122297 |

Magnetic stirring bar retrievers, PE



With a permanent magnet on one end and holding ring on the other one. Magnet is completely encapsulated in a PE mantle.

| Length mm | PU | Cat. No. |
|--------------|----|----------|
| 300 | 1 | 318293 |
| 450 | 1 | 318393 |

