



ADAPTORS FOR BOTTLES GAS WASHING, DRECHSLER, SJ 29/32

AUFSÄTZE FÜR GASWASCHFLASCHEN

nach Drechsler, NS 29/32

Cat. No. 8163

Catalog No. Katalog-Nr.	d (mm)	h1 (mm)	h2 (mm)	l (mm)	Packing (pcs) Packung (Stück)
1632441634080	8	180	50	40	10



BOTTLES GAS WASHING WITH SINTERED DISC

Typical use: cleaning (washing) gas.

FRITTEN – WASCHFLASCHE

Typische Verwendung: Reinigung
(Waschen) von Gas.

Cat. No. 815

Catalog No. Katalog-Nr.	Granulity Korngröße	Capacity (ml) Kapazität	h (mm)	h1 (mm)	SJ2	Packing (pcs) Packung (Stück)
1632442281102	s1	100	20	250	29/32	1
1632442281202	s2	100	20	250	29/32	1
1632442281103	s1	250	30	250	45/40	1
1632442281203	s2	250	30	250	45/40	1

