



BOWL WITH EDGE

Stainless steel 1.4301

Content in ml	200	400	1000	2000	3000	4000	5000	6000	8000	11000
Height in mm	65	45	90	115	120	120	135	135	160	180
Ø mm	105	155	175	217	245	286	315	320	350	395
	80.705.100	80.705.150	80.705.180	80.705.240	80.705.260	80.705.270	80.705.300	80.705.330	80.705.360	80.705.410



KIDNEY DISH

Stainless steel 1.4301

Length in mm	170	250	275
Content in ml	250	500	750
	80.708.170	80.708.250	80.708.275



EVAPORATING DISH

Stainless steel 1.4301

Content in ml	100	175	400
Height in mm	25	25	25
Ø in mm	75	100	150
	80.870.070	80.870.100	80.870.150



BEAKER

Molybdenum-alloyed, stainless steel 1.4404

Content in ml	100	250	500	1000	2000	3000	5000
Height in mm	65	80	110	130	180	180	210
Ø in mm	55	65	80	100	120	150	180
A flanged edge, smooth	80.075.100	80.075.250	80.075.500	80.075.001	80.075.002	80.075.003	80.075.005
B with spout	80.070.100	80.070.250	80.070.500	80.070.001	80.070.002	80.070.003	80.070.005
C smooth flanged edge with handle	80.078.100	80.078.250	80.078.500	80.078.001	80.078.002	80.078.003	80.078.005
D with sprout and handle	80.077.100	80.077.250	80.077.500	80.077.001	80.077.002	80.077.003	80.077.005