



Glass measuring cylinder, class B

Code **DBG525, DBG526, DBG527, DBG528**

- With spout
- Hexagonal base
- White graduation



Technical specifications

| Reference | Capacity (mL) | Tolerance (mL) | Graduation (mL) | Height (mm) | Diameter (mm) |
|-----------|---------------|----------------|-----------------|-------------|---------------|
| DBG525 | 15 | $\pm 0,2$ | 0,2 | 165 | 15 |
| DBG526 | 30 | $\pm 0,5$ | 0,5 | 210 | 20 |
| DBG527 | 60 | $\pm 1,0$ | 1,0 | 210 | 27 |
| DBG528 | 125 | $\pm 1,0$ | 1,0 | 268 | 32 |