

TAPERED STOPPERS - STANDARD DIN SIZES (according to DIN-12 871)

Verisilic grade:

Contains pure silicone. Taste- and odor-free. Amber colored.

Chemically inert, biologically neutral. Resistance to aging, ozone and repeated sterilization.

Temperature resistance from -60°C to 200°C (peaks acceptable up to 220°C).

Food-grade.



WITH STARTER HOLES

Type	DIN	Dimensions (mm)			SOLID				WITH STARTER HOLES			
		Dia. base	Dia. top	Height	Pkg. Qty.	Item no.	Other colors	Item no.	Number	Hole mm	Pkg. Qty.	Item no.
6D	6.5/15	3.5	6.5	15	500	263006						
8D		4	8	20	500	263008	Red	263007				
10D	9/20	5	9	20	500	263010	Blue	263011				
13D	12/20	8	12	20	500	263013	Grey	263014				
15D	14.5/20	10.5	14.5	20	500	263015						
16D	16.5/20	12.5	16.5	20	250	263016						
18D	18/20	14	18	20	250	263018			2	3	2500	263019
21D	22/25	17	22	25	100	263021			2	4	1400	263023
24D	24/30	18	24	30	50	263024						
27D	27/30	21	27	30	50	263027						
29D	29/30	23	29	30	50	263029						
31D	32/30	26	32	30	25	263031						
34D	35/30	29	35	30	25	263034						
38D	38/35	31	38	35	25	263038						
41D		34	41	35	25	263041						
45D	44/40	36	44	40	25	263045						
49D	49/40	41	49	40	20	263049						
54D	55/40	47	55	40	20	263054						
60D	59.5/45	50.5	59.5	45	20	263060						
65D	65/45	56	65	45	10	263065						
70D	70/50	60	70	50	5	263070						
75D	75.5/55	64.5	75.5	55	5	263075						
83D	83/60	71	83	60	5	263083						
92D	92/65	79	92	65	5	263092						
100D	100/65	87	100	65	5	263100						
107D	107/65	94	107	65	1	263107						

These items are not in stock and must be ordered specially.

Other colors can be supplied on request.