

SCREW FITTING

Screw fittings and fittings with GL thread for hard piping or pipes made of glass, metal, or plastic). For attachment of probes, thermometers, or cables in the reaction vessels. Ideal for the aggressive environments (aggressive vapors).


Green sleeve nut (PP) for the bottle GL45 thread and PP body. Maximum external diameter of the inserted pipe is 8,5mm. Limited chemical resistance, for operational temperatures up to +110°C. Suitable for inflow/outflow of liquids, leading-in the piping, pipes or probes into the vessels.

VERSCHRAUBUNG

Verschraubungen und Fittings mit GL-Gewinde für harte Rohrleitungen oder Rohre aus Glas, Metall oder Kunststoff. Befestigung von Sonden, Thermometern, Rohren oder Kabeln in Laborflaschen. Ideal für aggressive Umgebungen (mit aggressiven Dämpfen). Grüne Überwurfmutter aus PP für Flaschengewinde GL 45 und Körper aus PP. Höchstmöglicher Außendurchmesser des eingelegten Rohres 8,5mm.

Beschränkte chemische Beständigkeit, für Arbeitstemperaturen bis +110°C. Verwendung: Abführung oder Zuführung von Flüssigkeiten. Einführung der Rohrleitung, des Rohres oder der Sonden in die Behälter.

| Catalog No. Katalog-Nr. | Necks GL Hals-GL | Code | Packing (pcs) Packung (Stück) |
|----------------------------|---------------------|----------|----------------------------------|
| N555D612-08 | 2 x 14 | D 612-08 | 1 |
| N555D613-08 | 3 x 14 | D 613-08 | 1 |



SCREW CAP

with two hose connections

Green screw cap (PP) for the bottle GL45 thread plus body with hose connection (PP). The upper hose connections are bent, the lower ones are straight. Limited chemical resistance, for operational temperatures up to +110°C.

Used for inflow/outflow of liquids. Suitable for hoses (Vitron, Tygon, silicone).

SCHRAUBVERSCHLUSS - KAPPE

mit 2 Oliven für Schläuche

Grüner Kunststoffverschluss aus PP für Flaschengewinde GL 45 und Körper mit Oliven für Schläuche aus PP. Gebogene Oliven auf der Oberseite, gerade Oliven auf der Unterseite. Beschränkte chemische Beständigkeit, für Arbeitstemperaturen bis +110°C.

Verwendung: Abführung oder Zuführung von Flüssigkeiten. Für Schläuche (aus Vitron, Tygon, Silikon).

| Catalog No. Katalog-Nr. | Number of hose connectors Anzahl von Schläuchen | For tubing I.D. Für Rohr | Bore of hose connectors Innendurchmesser von Schläuchen | Code | Packing (pcs) Packung (Stück) |
|----------------------------|--|-----------------------------|--|----------|----------------------------------|
| N555D800-24 | 2 | 6-9 | 6 | D 800-24 | 1 |
| N555D800-36 | 3 | 6-9 | 6 | D 800-36 | 1 |
| N555D800-48 | 2 | 7-11 | 7 | D 800-48 | 1 |

