

Büchner, Guko adapters

Guko Adapters

Code PNJ003, PNJ004, PNJ005

- Made of natural rubber in red colour
- For Kitasato+Büchner connections



Code PNJ004

Code	Size	Upper external Ø	Lower external Ø	Height	Thickness
PNJ003	1	46 mm	28 mm	30,5 mm	3 mm
PNJ004	2	56 mm	38 mm	35 mm	3 mm
PNJ005	3	68 mm	47 mm	41 mm	4 mm

Set of 7 Guko adapters

Code PNJ006, PNJ007

- Made of rubber/silicone in blue colour
- For Kitasato+Büchner connections
- The dimensions of each adapter are:



N°	Ø ext larger/smaller	Ø int larger/smaller	Height
1	21x11 mm	17x7 mm	21 mm
2	27x16 mm	22x11 mm	22 mm
3	37x22 mm	31x16 mm	25 mm
4	46x29 mm	39x22 mm	29 mm
5	58x38 mm	50x30 mm	35 mm
6	69x45 mm	60x36 mm	40 mm
7	86x57 mm	75x46 mm	45 mm

Code	Description
PNJ006	Set of 7 adapters, rubber
PNJ007	Set of 7 adapters, silicone