



592



592D+595



593



592DCP



598CP

### ADJUSTABLE STOOLS

Seat Ø330, 60mm thick, high density foam 40kg/m<sup>3</sup>, flame retardant and washable upholstery (M1 class). Foot operated height adjustment by gaz spring from 600 to 700mm. Five legs base Ø535mm, made in chromed steel, 5 castors Ø50mm with non marking type. Round tray Ø350mm for height adjustment by foot.

Seat color : ref.591CP : Black – ref. 591CPMA : Mandarine – ref.

591CPMI : Mirabelle – ref. 591CPC : Blue Cobalt – ref. 591CPE : Blue Egee.

### ADJUSTABLES STOOLS

Seat vinyl without seam Ø 330 mm. Anti-fire quality M2. Gaz spring handfull adjustment (5 years warranted covered). 5 twin-castors 50 mm diameter or plastic skids.

#### 5 legs aluminium base

| REF.    | Seat  | Height       | Base   |
|---------|---|--------------|--------|
| 592     | 330 mm  | 520 à 700 mm | 520 mm |
| 592CP   | foot control ht. adjustable                             |              |        |
| 592D    | with backrest   |              |        |
| 592DCP  | with backrest, foot control                             |              |        |
| 592P    | 5 plastic skids   |              |        |
| 592PCP  | 5 plastic skids, foot control ht adjustable             |              |        |
| 592DP   | with backrest + 5 plastics skids                        |              |        |
| 592DPCP | backrest + 5 plastics skids, foot control ht adjustable |              |        |

#### 5 legs nylon base

| REF.    | Seat   | Height       | Base   |
|---------|--|--------------|--------|
| 593     | 330 mm   | 520 à 700 mm | 520 mm |
| 593CP   | foot control ht. adjustable                                  |              |        |
| 593D    | with backrest  |              |        |
| 593DCP  | with backrest  |              |        |
| 593P    | 5 plastic skids  |              |        |
| 593PCP  | 5 plastic skids, foot control ht adjustable                  |              |        |
| 593DP   | with backrest + 5 plastics skids                             |              |        |
| 593DPCP | with backrest + 5 plastics skids, foot control ht adjustable |              |        |

#### Adjustable foot rest for stool (not for foot adjustment)

595 adjustable in height

#### 5 legs stainless steel base. Hydraulic foot adjustment

| REF.    | Diam.Seat                | Height       | Base   |
|---------|--------------------------|--------------|--------|
| 598CP   | 330 mm skids             | 500 à 700 mm | 520 mm |
| 598DCP  | with backrest            |              |        |
| 598CPR  | 330 mm castors           | 520 à 720 mm | 520 mm |
| 598DCPR | with backrest 0 à 720 mm |              |        |
| 598DCPR | Idem avec dossier        |              |        |