



CASSETTES AND ACCESSORIES

Code **BPJ001_011**, **BPJ016_020**, **BPJ040_044**

Cassettes for embedding in paraffin, made of solvent-resistant plastic, compatible with major processing systems. Several models available: universal, for standard samples, and for biopsy, with finer perforations for small tissues. They come with an integrated lid or without a lid, and in various colours.

Additionally, accessories and complementary items are available to adapt the use of the cassettes to different laboratory needs.

Universal cassettes with lid - Code **BPL001_005**



- Made of polyoxymethylene (POM)
- Fully resistant to organic solvents used in histological processing
- Feature a large 45° angled labelling surface for easy manual or automated identification
- 1x5 mm perforations that allow optimal reagent circulation without compromising sample integrity

Code	Colour	Approx. Dimensions
BPJ001	White	40x28x6 mm
BPJ002	Blue*	
BPJ003	Pink*	
BPJ005	Green*	

*While stock last


Universal cassettes without lid - Code **BPJ011**



- Made of polyoxymethylene (POM)
- Fully resistant to organic solvents used in histological processing
- Feature a large 45° angled labelling surface for easy manual or automated identification
- 2 mm diameter perforations that allow optimal reagent circulation without compromising sample integrity

Code	Colour	Approx. Dimensions
BPJ011	White	40x28x6 mm

Biopsy cassettes with lid - Code **BPJ006_010**



- Made of polyoxymethylene (POM)
- Fully resistant to organic solvents used in histological processing
- Feature a large 45° angled labelling surface for easy manual or automated identification
- Microperforations of 1x1 mm that allow uniform reagent flow while maintaining the integrity of very small samples

Code	Colour	Approx. Dimensions
BPJ006	White	40x28x6 mm
BPJ008	Pink*	
BPJ009	Yellow*	
BPJ010	Green*	

*While stock last



CASSETTES AND ACCESORIES

Code BPJ001_011, BPJ016_020, BPJ040_044

Sponges for cassettes - Code BPJ044



- Foam pads designed to secure small samples during histological processing
- Keep biopsies in place within the cassette, preventing the loss of fragments during processing
- Made of cellular polyurethane polyester foam
- Contain microcavities that allow uniform reagent flow around the tissue, without compromising the integrity of the sample

Code	Colours	Approx. Dimensions
BPJ044	White, blue	28x25x2 mm

Embedding rings - Code BPJ016_020



- Made of high-quality ABS plastic
- Feature a large labelling surface
- Reinforced with two bars in the central hole to prevent bending or deformation

Code	Colour	Ring dimensions	Base dimensions
BPJ016	White	25x25 mm	40x28 mm
BPJ017	Yellow*		
BPJ018	Green*		
BPJ019	Pink*		
BPJ020	Blue*		

*While stock last

Embedding molds - Code BPJ040_043



- Made of stainless steel
- Available in various models depending on the size of the resulting block, to accommodate different types of samples

Code	Block dimensions
BPJ040	7x7x6 mm
BPJ041	15x15x5 mm
BPJ042	24x24x6 mm
BPJ043	30x24x6 mm