



## VOLUMETRIC FLASK WITH SCREW CAP, PP

Code DDP007 to DDP012

- Made of high-quality transparent polypropylene (PP)
- Narrow neck design for precise pouring
- Autoclavable



### Technical specifications

Code	Capacity (mL)	Outer Ø base (mm)	Height (mm)	GL (Thread)
DDP007	25	36	121	18
DDP008	50	46	143	18
DDP009	100	58	174	18
DDP010	250	82	225	25
DDP011	500	100	265	25
DDP012	1000	120	325	25