



Flotation bath for tissue sections

CODE ZFP008

Thermostatic bath for paraffin sections. Equipment with new design and materials, high functional stability and easy to use, it is used in Histology and Pathology for tissues treatment.



- Wide panel control, which is easy to understand and use.
- At the appropriate temperature the thermostatic chamber permits the separation and subsequent application on the microscope slide of the strip obtained with the microtome.
- Black surface of the chamber for a good observation of the section.
- Wide thermal limit to get an excellent drying of the preparation.

Technical specifications

Code	ZFP008
Capacity	2L (approx.)
Temperature range	Room temperature to 90°C
Precision	-3 / +5 °C
Dimensions (chamber)	240x180x54 mm
Dimensions (equipment)	340x350x105 mm
Supply	AC 220±10% (50Hz)
Consume	550 w
Fuses	6A
Weight	4.5 kg