

Chromatography Column With Sintered/Fritted Disc

- With Socket
- Fritted Disc Porosity 0
- With PTFE Key Stopcock, 2mm Bore
- Made from Boro 3.3 Glass

Item Code	I.D. (mm)	Length (mm)	Capacity (ml)	Stopcock	N/S	Pack of
221.212.01	10	100	8	PTFE	14/23	1
221.212.02	10	200	15	PTFE	14/23	1
221.212.03	15	200	35	PTFE	14/23	1
221.212.04	10	300	23	PTFE	14/23	1
221.212.05	20	400	125	PTFE	29/32	1
221.212.07	30	600	430	PTFE	29/32	1

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

PTFE

1 x
SOCKET



ENQUIRE



Chromatography Column With Sintered/Fritted Disc

- With Socket and Cone
- With PTFE Key Stopcock and Suction Tube (GL 14)
- Fritted Disc Porosity 0
- Made from Boro 3.3 Glass

Item Code	I.D. (mm)	Length (mm)	Capacity (ml)	Socket & Cone	Pack of
212.202.02	15	200	35	14/23	1
212.202.04	20	400	120	14/23	1
212.202.04D	30	400	280	29/32	1
212.202.06D	30	600	430	29/32	1

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

PTFE

1 x
SOCKET



ENQUIRE

