

## Chromatography Column With Sintered/Fritted Disc

- With Socket
- Fritted Disc Porosity 0
- With PTFE Key Stopcock, 2mm Bore
- Made from Boro 3.3 Glass

| Item Code  | I.D. (mm) | Length (mm) | Capacity (ml) | Stopcock | N/S   | Pack of |
|------------|-----------|-------------|---------------|----------|-------|---------|
| 221.212.01 | 10        | 100         | 8             | PTFE     | 14/23 | 1       |
| 221.212.02 | 10        | 200         | 15            | PTFE     | 14/23 | 1       |
| 221.212.03 | 15        | 200         | 35            | PTFE     | 14/23 | 1       |
| 221.212.04 | 10        | 300         | 23            | PTFE     | 14/23 | 1       |
| 221.212.05 | 20        | 400         | 125           | PTFE     | 29/32 | 1       |
| 221.212.07 | 30        | 600         | 430           | PTFE     | 29/32 | 1       |

**3.3**  
BORO  
GLASS

WHITE  
PRINTING

Chemical  
Resistant

PTFE

1 x  
SOCKET



ENQUIRE



## Chromatography Column With Sintered/Fritted Disc

- With Socket and Cone
- With PTFE Key Stopcock and Suction Tube (GL 14)
- Fritted Disc Porosity 0
- Made from Boro 3.3 Glass

| Item Code   | I.D. (mm) | Length (mm) | Capacity (ml) | Socket & Cone | Pack of |
|-------------|-----------|-------------|---------------|---------------|---------|
| 212.202.02  | 15        | 200         | 35            | 14/23         | 1       |
| 212.202.04  | 20        | 400         | 120           | 14/23         | 1       |
| 212.202.04D | 30        | 400         | 280           | 29/32         | 1       |
| 212.202.06D | 30        | 600         | 430           | 29/32         | 1       |

**3.3**  
BORO  
GLASS

WHITE  
PRINTING

Chemical  
Resistant

PTFE

1 x  
SOCKET



ENQUIRE

