

# Funnels



## Funnels, Pressure Equalizing, Pear Shape

- Made from DIN ISO 3585, BORO 3.3 GLASS
- With Stem inside cone

Item Code	Capacity (ml)	N/S	Bore Size	Stopcock Type	Pack of
170.202.01	50	14/23	2.5	Solid Glass	1
170.202.02A	100	14/23	2.5	Solid Glass	1
170.202.04A	250	29/32	4	Solid Glass	1
170.202.05A	500	29/32	4	Solid Glass	1

3.3  
BORO  
GLASS

Chemical  
Resistant

WHITE  
PRINTING

GLASS

PE

WITH CONE  
& SOCKET



ENQUIRE



## Funnel, Pressure Equalizing, Cylindrical

- Made from DIN ISO 3585, BORO 3.3 GLASS
- With Stem inside cone

Item Code	Capacity (ml)	N/S	Bore (mm)	Stopcock Type	Pack of
168.202.01	50	14/23	2.5	Solid Glass	1
168.202.01A	50	29/32	2.5	Solid Glass	1
168.202.02A	100	14/23	2.5	Solid Glass	1
168.202.03A	100	29/32	2.5	Solid Glass	1
168.202.04A	250	29/32	4	Solid Glass	1
168.202.05A	500	29/32	4	Solid Glass	1

3.3  
BORO  
GLASS

Chemical  
Resistant

WHITE  
PRINTING

GLASS

PE

WITH CONE  
& SOCKET



ENQUIRE