


Funnels



Funnel

- Made from ASTM E-438 Type 1, CLASS A, Boro 3.3 glass
- Complies with DIN ISO 4798 & USP standards
- With 60° bowl angle for precise fitting of filter paper

Item Code	Diameter (mm)	Pack of
238.202.01*	25	10
238.202.02	35	10
238.202.03*	38	10
238.202.04*	40	10
238.202.05*	50	10
238.202.06	55	10
238.202.07*	64	10
238.202.08*	70	10
238.202.09	75	10
238.202.10*	80	10
238.202.11	100	10
238.202.12	125	10
238.202.13	150	10

3.3
BORO
GLASS



* Non DIN



ENQUIRE



Funnel, Long stem

- With 60° bowl angle for precise fitting of filter paper
- Manufactured from chemical and heat resistant borosilicate glass
- Stem length approximately 150mm

Item Code	Dia (ml)	Stem O.D. (mm)	To suit filler paper dia (mm)	Pack of
238.215.08	70	8	70-90	6
238.215.10	80	10	90-150	6
238.215.11A	105	10	150-185	6

3.3
BORO
GLASS



ENQUIRE