

Cylinder-PP, Measuring, Hexagonal

• Hexagonal base provides excellent stability to the cylinder

Item Code	Capacity (ml)	Least Count (ml)	Pack of
177.303.01	10	0.2	24
177.303.02	25	0.5	12
177.303.03	50	1	12
177.303.04	100	1	12
177.303.05	250	2	6
177.303.06	500	5	4
177.303.07	1000	10	2



ENQUIRE



Cylinder, Tapped Density Cylinder

• Complies USP 616

• Used in tapped density apparatus to measure bulk density of powders

Item Code	Capacity (ml)	Pack of
139.520.01A	100 & 250	1

Cylinder, Tapped Density Cylinder

Item Code	Capacity (ml)	Pack of
139.521.01A	100	2

Cylinder, Tapped Density Cylinder

Item Code	Capacity (ml)	Pack of
139.521.05A	250	2

3.3
BORO
GLASS



ENQUIRE



ENQUIRE



ENQUIRE



Cylinder, Brush

Item Code	Type	Material of Bristle	Pack of
593.303.05	Cylinder Brush	Nylon	50

Nylon



ENQUIRE

