

Joints



Cone, Unprinted

Item Code	Socket Size	Approx. O.D. of tube (mm)	Total Length (mm)	Pack of
002.202.02	10/19	8	120±10	10
002.202.03	12/21	9	120±10	10
002.202.04	14/23	13	120±10	10
002.202.05	19/26	16	125±10	10
002.202.06	24/29	22	135±10	10
002.202.07	29/32	26	135±10	10
002.202.08	34/35	30	150±10	10
002.202.09*	40/38	36	150±10	10
002.202.10	45/40	40	150±10	10
002.202.11*	50/42	45	150±10	10
002.202.12*	55/44	50	150±10	10

3.3
BORO
GLASS



* Non DIN



ENQUIRE

Cone, With drip tip, Unprinted

Item Code	Socket Size	Approx. O.D. of tube (mm)	Total Length (mm)	Pack of
003.202.03	12/21	9	120±10	10
003.202.04	14/23	13	120±10	10
003.202.05	19/26	16	120±10	10
003.202.06	24/29	22	150±15	10
003.202.07	29/32	26	150±15	10
003.202.08	34/35	30	150±15	10
003.202.09*	40/38	36	150±15	10
003.202.10	45/40	40	150±15	10
003.202.11*	50/42	45	150±15	10
003.202.12*	55/44	50	150±15	10

3.3
BORO
GLASS



* Non DIN



ENQUIRE

Socket, Double

• Complies with DIN 12249 standard

Item Code	Socket Size	Approx. O.D. of tube (mm)	Total Length (mm)	Pack of
006.202.02*	10/19	14	135±2	10
006.202.03	14/23	18	135±2	10
006.202.04	19/26	22	140±2	10
006.202.05	24/29	28	170±5	10
006.202.06	29/32	32	170±5	10
006.202.07	34/35	38	180±5	10

3.3
BORO
GLASS



* Non DIN



ENQUIRE