

Volumetric Flasks, Class A

- Complies with DIN ISO 1042 Standards
- With Batch certificate
- With PE Stopper

Item Code	Capacity (ml)	N/S	Tolerance (±)	Pack of
130.202.01	5	10/19	0.040	2
130.202.02	10	10/19	0.040	2
130.202.02A	20	10/19	0.040	2
130.202.03	25	10/19	0.040	2
130.202.04	50	12/21	0.060	2
130.202.04B	50	14/23	0.100	2
130.202.05B	100	12/21	0.100	2
130.202.05	100	14/23	0.100	2
130.202.06	200	14/23	0.150	2
130.202.07	250	14/23	0.150	2
130.202.08	500	19/26	0.250	2
130.202.09	1000	24/29	0.400	2
130.202.10	2000	29/32	0.600	1
130.202.11	5000	34/35	1.200	1



ENQUIRE



3.3
BORO
GLASS

Chemical
Resistant

BLUE
PRINTING

A
CLASS

PE

BLUE
RING

BATCH
CERTIFICATE

Trapezoidal Volumetric Flask, CLASS A, Clear Glass *New*

- Boro 3.3 glass
- With Batch Certificate
- With PE Stopper

Item Code	Capacity (ml)	N/S	Pack of
135.202.00	1	7/16	2
135.202.00B	2	7/16	2
135.202.01	5	10/19	2
135.202.02	10	10/19	2
135.202.03	20	10/19	2
135.202.04	25	10/19	2
135.202.05	50	12/21	2



ENQUIRE



3.3
BORO
GLASS

Chemical
Resistant

BLUE
PRINTING

A
CLASS

PE

BLUE
RING

BATCH
CERTIFICATE