


Adapter



Adapter, Vacuum Receiver, Angled, 105°

- Complies with DIN 12594 standard
- With drip tip and delivery stem inside the cone
- With a screw thread and bend with Vacuum Connection
- Used for distribution of liquid from the condenser into collection vessel

Item Code	Socket Size	Cone Size	Pack of
022.201.01	14/23	14/23	1
022.201.02*	14/23	19/26	1
022.201.03	19/26	19/26	1
022.201.04*	19/26	24/29	1
022.201.05*	24/29	24/29	1
022.201.06*	24/29	29/32	1
022.201.07	29/32	29/32	1

*Non DIN/ISO

ASTM Sizes

Item Code	Cone Size to fit flask	Cone Size to fit Condenser	Pack of
022.471.01S	14/20	14/20	1
022.471.03S	19/22	19/22	1
022.471.05S	24/40	24/40	1

3.3
BORO
GLASS



1 x
SOCKET

1 x
CONE



ENQUIRE



ENQUIRE

Adapter, Receiver With Side Socket

- Socket inclined at 105° angle
- With drip tip inside the cone

Item Code	Cone Size to fit flask	Cone Size to fit Condenser	Pack of
023.202.01	14/23	14/23	1
023.202.02	19/26	19/26	1
023.202.03	19/26	24/29	1
023.202.04	24/29	24/29	1
023.202.05	29/32	29/32	1

3.3
BORO
GLASS



2 x
SOCKET

1 x
CONE



ENQUIRE



Adapter, Distilling Cow Receiver

- Four cones encased with drip tip and one socket joint
- Flask 100 ml

Item Code	Socket Size	Cone Size	Pack of
024.475.01	14/20	14/20	1

3.3
BORO
GLASS



1 x
SOCKET

4 x
CONE



ENQUIRE

