



Cone, Unprinted

| Item Code | Socket Size | Approx. O.D. of tube (mm) | Total Length (mm) | Pack of |
|-------------|-------------|---------------------------|-------------------|---------|
| 002.202.02 | 10/19 | 8 | 120±10 | 10 |
| 002.202.03 | 12/21 | 9 | 120±10 | 10 |
| 002.202.04 | 14/23 | 13 | 120±10 | 10 |
| 002.202.05 | 19/26 | 16 | 125±10 | 10 |
| 002.202.06 | 24/29 | 22 | 135±10 | 10 |
| 002.202.07 | 29/32 | 26 | 135±10 | 10 |
| 002.202.08 | 34/35 | 30 | 150±10 | 10 |
| 002.202.09* | 40/38 | 36 | 150±10 | 10 |
| 002.202.10 | 45/40 | 40 | 150±10 | 10 |
| 002.202.11* | 50/42 | 45 | 150±10 | 10 |
| 002.202.12* | 55/44 | 50 | 150±10 | 10 |

3.3
BORO
GLASS



* Non DIN



ENQUIRE



Cone, With drip tip, Unprinted

| Item Code | Socket Size | Approx. O.D. of tube (mm) | Total Length (mm) | Pack of |
|-------------|-------------|---------------------------|-------------------|---------|
| 003.202.03 | 12/21 | 9 | 120±10 | 10 |
| 003.202.04 | 14/23 | 13 | 120±10 | 10 |
| 003.202.05 | 19/26 | 16 | 120±10 | 10 |
| 003.202.06 | 24/29 | 22 | 150±15 | 10 |
| 003.202.07 | 29/32 | 26 | 150±15 | 10 |
| 003.202.08 | 34/35 | 30 | 150±15 | 10 |
| 003.202.09* | 40/38 | 36 | 150±15 | 10 |
| 003.202.10 | 45/40 | 40 | 150±15 | 10 |
| 003.202.11* | 50/42 | 45 | 150±15 | 10 |
| 003.202.12* | 55/44 | 50 | 150±15 | 10 |

3.3
BORO
GLASS



* Non DIN



ENQUIRE



Socket, Double

• Complies with DIN 12249 standard

| Item Code | Socket Size | Approx. O.D. of tube (mm) | Total Length (mm) | Pack of |
|-------------|-------------|---------------------------|-------------------|---------|
| 006.202.02* | 10/19 | 14 | 135±2 | 10 |
| 006.202.03 | 14/23 | 18 | 135±2 | 10 |
| 006.202.04 | 19/26 | 22 | 140±2 | 10 |
| 006.202.05 | 24/29 | 28 | 170±5 | 10 |
| 006.202.06 | 29/32 | 32 | 170±5 | 10 |
| 006.202.07 | 34/35 | 38 | 180±5 | 10 |

3.3
BORO
GLASS



* Non DIN



ENQUIRE