

Condensers Allihin for Soxhlet Apparatus

- Manufactured from ASTM E-438, Type-1, Class-A boro 3.3 Glass
- Condenser with glass screw thread connector

Item Code	Cone	Effective Length (mm)	Extractor (ml)	Connector	Pack of
191.202.01	45/40	300	100	Glass	1
191.202.01A	45/40	400	250	Glass	1
191.202.02	60/46	400	500	Glass	1
191.202.03	71/55	400	1000	Glass	1

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant



ENQUIRE



Simple Distillation Equipment Set

- Round Bottom Flask 250ml, NS 29/32
- Adapter, Still head plain socket 29/32, Cone 29/32, 29/32
- Liebig Condenser 250mm 29/32 with Glass hose connector
- Adapter Receiver 29/32, short stem
- Adaptor Thermometer with compression cap 29/32
- Thermometer 10° to 110° C

Item Code	Cone	Effective Length (mm)	Extractor (ml)	Connector	Pack of
191.181.01	45/40	300	100	Screw Thread	1
191.181.01A	45/40	400	250	Screw Thread	1
191.181.02	60/46	400	500	Screw Thread	1
191.181.03	71/55	400	1000	Screw Thread	1

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

SCREW
CAP



ENQUIRE

