

Burette



Automatic Burette, Unserialized

- With Batch Certificate
- Compliance with ASTM E-287 Standard
- With Intermediate PTFE Key Stopcock And With Needle Valve Stopcock

Item Code	Capacity (ml)	Division	Tolerance ($\pm$ ml)	Pack of
118.456.01	10	0.05	0.02	1
118.456.02	25	0.10	0.03	1
118.456.03	50	0.10	0.05	1

3.3
BORO
GLASS

BLUE
PRINTING

Chemical
Resistant

PTFE

NEEDLE

AS
CLASS

BATCH
CERTIFICATE



ENQUIRE

Accessories like bottle & rubber below order separately (For reference go to page no 77)
Ps: These burettes can also be supplied in amber color. Please ask for the same separately



Schilling Burette with Schellbach Stripe, Class B

- Sub. division & tolerance as per DIN/ISO 385 Standard
- With Plastic Bottle 1000ml
- With Plastic Clamp, Stand and PTFE Pipe
- Schellbach Stripe for Precise Liquid Dispensing and Accurate Meniscus Reading

Item Code	Capacity (ml)	Division	Tolerance ($\pm$ ml)	Pack of
118.224.01	10	0.05	0.05	1
118.224.02	25	0.10	0.10	1
118.224.03	50	0.10	0.10	1

Schilling Burette with Schellbach Stripe, Amber

Item Code	Capacity (ml)	Division	Tolerance ($\pm$ ml)	Pack of
118.237.01AMB	10	0.05	0.05	1
118.224.02AMB	25	0.10	0.10	1
118.224.03AMB	50	0.10	0.10	1



Schellbach
Stripe

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

PTFE

B
CLASS



ENQUIRE



ENQUIRE