

Flask



Flask, Boiling, Round Bottom with Narrow Neck

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- Complies with ISO 1773 standard
- With beaded rim


Item Code	Capacity (ml)	Body Dia (mm)	Neck Dia (mm)	Height (mm)	Pack of
233.202.01	50	51	26	95	10
233.202.02	100	64	26	110	10
233.203.02*	150	74	28	118	10
233.202.03	250	85	34	143	10
233.202.04	500	105	34	168	10
233.202.05	1000	131	42	200	10
233.202.06	2000	166	50	240	6
233.202.07*	3000	185	50	260	1
233.202.08*	5000	223	50	305	1
233.202.09	10000	279	65	380	1
233.202.10*	20000	350	75	510	1

3.3
BORO
GLASS



ENQUIRE

*Non DIN



Flask, Boiling, Round Bottom with Wide Neck

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- Complies with ISO 24450 standard
- With beaded rim

Item Code	Capacity (ml)	Body Dia (mm)	Neck Dia (mm)	Height (mm)	Pack of
234.202.01	50	51	34	95	10
234.202.02	100	64	34	110	10
234.202.03	250	85	50	143	10
234.202.04	500	105	50	168	10
234.202.05	1000	131	50	200	10
234.202.07	2000	166	76	240	6

3.3
BORO
GLASS



ENQUIRE