

# Cylinder

GLASSCO measuring cylinders are made from ASTM E – 438, TYPE-1, CLASS-A BORO 3.3 glass that meet ISO 4788 standard. The calibration is done on automatic computer controlled machines and are calibrated to contain (TC, In) according to the ISO & ASTM tolerance standards. All ISO Cylinders come with **[DE-M]** Mark printed on it.



## Measuring Cylinder with Hexagonal Base, Class-A

- Certificate as per DIN EN ISO 4788
- With individual work certificate

Item Code	Capacity (ml)	Sub Division (ml)	Tolerance ( $\pm$ ml)	Pack of
139.223.00A	5	0.1	0.05	1
139.223.01A	10	0.2	0.10	1
139.223.02A	25	0.5	0.25	1
139.223.03A	50	1.0	0.50	1
139.223.04A	100	1.0	0.50	1
139.223.05A	250	2.0	1.00	1
139.223.06A	500	5.0	2.50	1
139.223.07A	1000	10.0	5.00	1
139.223.08A	2000	20.0	10.00	1

3.3  
BORO  
GLASS

BLUE  
PRINTING

Chemical  
Resistant

A  
CLASS

HEXA  
BASE

INDIVIDUAL  
CERTIFICATE



ENQUIRE



## Measuring Cylinder with Hexagonal Base, Class-A

- Certificate as per DIN EN ISO 4788
- With conformity Batch certificate

Item Code	Capacity (ml)	Sub Division (ml)	Tolerance ( $\pm$ ml)	Pack of
139.202.00A	5	0.1	0.05	2
139.202.01A	10	0.2	0.10	2
139.202.02A	25	0.5	0.25	2
139.202.03A	50	1.0	0.50	2
139.202.04A	100	1.0	0.50	2
139.202.05A	250	2.0	1.00	2
139.202.06A	500	5.0	2.50	2
139.202.07A	1000	10.0	5.00	1
139.202.08A	2000	20.0	10.00	1

3.3  
BORO  
GLASS

BLUE  
PRINTING

Chemical  
Resistant

A  
CLASS

HEXA  
BASE

BATCH  
CERTIFICATE



ENQUIRE

PS: To order B class low form cylinder go to page :