

Funnels



Funnel, Buchner, Euro Design

- Made of Porcelain
- Temperature upto 1000°C

Item Code	Dia (mm) OD x H	Pack of
527.303.10	20 x 64	10
527.303.11	51 x 89	10
527.303.12	67 x 110	6
527.303.13	84 x 143	4
527.303.14	105 x 160	2
527.303.15	125 x 195	2
527.303.16	142 x 202	1
527.303.17	170 x 221	1
527.303.18	210 x 270	1



ENQUIRE



Funnel, Buchner, With Sintered/Fritted Disc

- Made from ASTM E 438 type-1, Class A Boro 3.3 glass
- With fire polished rim and stem

Item Code	Capacity (ml)	Disk Dia (mm)	Porosity	Pack of
256.G00.01	35	30	0	2
256.G01.01	35	30	1	2
256.G02.01	35	30	2	2
256.G03.01	35	30	3	2
256.G04.01	35	30	4	2
256.G05.01	35	30	5	2
256.G00.02	80	40	0	2
256.G01.02	80	40	1	2
256.G02.02	80	40	2	2
256.G03.02	80	40	3	2
256.G04.02	80	40	4	2
256.G05.02	80	40	5	2
256.G00.03	200	65	0	2
256.G01.03	200	65	1	2
256.G02.03	200	65	2	2
256.G03.03	200	65	3	2
256.G04.03	200	65	4	2
256.G05.03	200	65	5	2
256.G00.04	500	90	0	2
256.G01.04	500	90	1	2
256.G02.04	500	90	2	2
256.G03.04	500	90	3	2
256.G04.04	500	90	4	2
256.G05.04	500	90	5	2
256.G00.05	1000	120	0	1
256.G01.05	1000	120	1	1
256.G02.05	1000	120	2	1
256.G03.05	1000	120	3	1
256.G04.05	1000	120	4	1
256.G05.05	1000	120	5	1

3.3
BORO
GLASS



ENQUIRE