

Flask, Filter

- With Glass Connector
- Made from heavy wall which can withstand high vacuum
- Printed in durable white colour enamel without Graduation

Item Code	Capacity (ml)	Pack of
074.202.00	50	2
074.202.01	100	2
074.202.02	250	2
074.202.03	500	2
074.202.04	1000	2
074.202.05	2000	1
074.202.08	5000	1

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

HEAVY
WALL



ENQUIRE

* Non DIN



Flask, Buchner, plain side arm and ground socket

- Complies with DIN ISO 3585 standard
- With precision-ground quick fit socket
- With traditional plain glass side-arm to accept 9mm bore vacuum tubing
- Heavy wall and base construction for vacuum work

Item Code	Capacity (ml)	Socket Size	Height (mm)	Approx Dia (mm)	Pack of
074.205.01B	100	19/26	107	70	2
074.205.01C	100	24/29	107	70	2
074.205.02	250	24/29	158	85	2
074.205.02A	250	29/32	158	85	2
074.205.03	500	24/29	188	105	2
074.205.03A	500	29/32	188	105	2
074.205.04	1000	24/29	238	138	2
074.205.04A	2000	29/32	238	138	2

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant

HEAVY
WALL



ENQUIRE

*Without Stopper

