

Adapter



Adapter, Socket to Cone With Side-Arm Screw Thread

- With Drip Tip
- Detachable plastic screw thread that makes it compatible for a rubber hose connection

Item Code	Socket Size	Cone Size	Pack of
039.202.01	14/23	14/23	1
039.202.02	14/23	19/26	1
039.202.03	14/23	24/29	1
039.202.04	14/23	29/32	1
039.202.05	19/26	19/26	1
039.202.06	19/26	24/29	1
039.202.07	19/26	29/32	1
039.202.08	24/29	24/29	1
039.202.09	24/29	29/32	1
039.202.10	29/32	29/32	1

3.3
BORO
GLASS



1 x
CONE

1 x
SOCKET

SCREW



ENQUIRE



Adapter, Vacuum or Gas 90°

- With connection for rubber hose

Item Code	Cone Size	Pack of
341.470.02	14/20	1
341.470.03	19/22	1
341.470.05A	24/40	1
341.470.05	29/42	1
341.470.04	35/20	1

3.3
BORO
GLASS



1 x
CONE



ENQUIRE



Adapter, Vacuum or Gas 90°

- With 2 connections for rubber hose

Item Code	Cone Size	Pack of
046.202.04	24/40	1

3.3
BORO
GLASS



1 x
CONE



ENQUIRE