

Flask



Flask, Kjeldahl

- Complies with DIN 12360 standard
- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- Suitable for determination of Nitrogen

Item Code	Capacity (ml)	D.(mm)	H.(mm)	Pack of
077.202.01	100	60	200	2
077.202.02*	300	85	295	2
077.202.03	500	101	300	2
077.202.04*	800	117	340	2

3.3
BORO
GLASS



* Non DIN



ENQUIRE



Flask, Kjeldahl

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- Used for the Determination of Nitrogen

Item Code	Capacity (ml)	Socket Size	Approx. Overall Height (mm)	Pack of
070.202.01	50	19/26	200	2
070.202.02	100	19/26	200	2
070.202.03	100	24/29	200	2
070.202.04*	300	24/29	295	2
070.202.05	500	24/29	300	2
070.202.06*	800	24/29	340	2

3.3
BORO
GLASS



* Non DIN



ENQUIRE



Flask, Brush

- Made of Nylon Bristles

Item Code	Type	Material of Bristle	Pack of
593.303.04	Flask Brush	Nylon	50



ENQUIRE