


JOINTS

- In compliance with ISO 383/DIN 12249 Standards.
- Unprinted sockets to avoid getting printed charred
- Smooth grinding surface finish to avoid leakage

Conical Joints

Size	Actual Dia of Wide End (mm)	Actual Dia of Narrow End (mm)	Nominal Length of Engagement (mm)
10/19	10.0	8.1	19
12/21	12.5	10.4	21
14/23	14.5	12.2	23
19/26	18.8	16.2	26
24/29	24.0	21.1	29
29/32	29.2	26.0	32
34/35	34.5	31.0	35
40/38*	40.0	36.2	38
45/40	45.0	41.0	40
50/42*	50.0	45.8	42
55/44*	55.0	50.6	44



* Non DIN

Socket, Unprinted

Item Code	Socket Size	Approx. O.D. of tube (mm)	Total Length (mm)	Pack of
001.202.02	10/19	13	120±10	10
001.202.03	12/21	16	120±10	10
001.202.04	14/23	18	120±10	10
001.202.05	19/26	22	125±10	10
001.202.06	24/29	28	135±10	10
001.202.07	29/32	32	135±10	10
001.202.08	34/35	38	150±10	10
001.202.09*	40/38	45	150±10	10
001.202.10	45/40	50	150±10	10
001.202.11*	50/42	55	150±10	10
001.202.12*	55/44	60	150±10	10



* Non DIN



ENQUIRE

