


Flask



Flask, Flat Bottom With Narrow Neck

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- With beaded rim
- Complies with DIN/ISO 1773

Item Code	Capacity (ml)	Body Dia (mm)	Neck Dia (mm)	Height (mm)	Pack of
235.202.01	50	51	26	90	10
235.202.02	100	64	26	106	10
235.204.02*	150	74	28	115	10
235.202.03	250	85	34	138	10
235.204.03*	300	93	34	143	10
235.202.04	500	105	34	163	10
235.202.05	1000	131	42	190	10
235.202.06	2000	166	50	230	6
235.202.07*	3000	185	50	250	1
235.202.08*	5000	233	50	290	1
235.202.09	10000	280	65	360	1
235.202.10*	20000	350	75	495	1

3.3
BORO
GLASS


WHITE
PRINTING

Chemical
Resistant



ENQUIRE

*Non DIN



Flask, Flat Bottom with Wide Neck

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- With beaded rim
- Complies with DIN/ISO 24450

Item Code	Capacity (ml)	Body Dia (mm)	Neck Dia (mm)	Height (mm)	Pack of
236.202.01	50	51	34	90	10
236.202.02	100	64	34	105	10
236.202.03	250	85	50	138	10
236.202.04	500	105	50	163	10
236.202.05	1000	131	50	190	10
236.202.06	2000	166	76	230	6

3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant



ENQUIRE