

Adapter, Still Head Plain

- Complies with DIN 12594 standard
- with drip tip, 75°
- with thermometer socket
- used for connecting flask to a condenser in a distillation assembly

Item Code	Socket Size	Cone Size to fit flask	Cone Size to fit Condenser	Pack of
028.202.01	14/23	14/23	14/23	1
028.202.02	14/23	19/26	19/26	1
028.202.03*	14/23	24/29	19/26	1
028.202.04	14/23	29/32	19/26	1
028.202.05*	14/23	34/35	19/26	1
028.202.06*	14/23	24/29	24/29	1
028.202.07*	14/23	34/35	24/29	1
028.202.08*	14/23	29/32	24/29	1
028.202.09*	14/23	29/32	29/32	1
028.202.10*	14/23	34/35	29/32	1

*Non DIN/ISO

ASTM Sizes

Item Code	Socket Size	Cone Size to fit flask	Cone Size to fit Condenser	Pack of
028.470.02	19/22	19/22	19/22	1
028.470.03	24/40	24/40	24/40	1
028.470.04	14/20	14/20	14/20	1

3.3
BORO
GLASS



1 x
SOCKET

2 x
CONE



ENQUIRE



ENQUIRE



Adapter, Splash Heads, Straight

- Pear shaped
- With drip tip
- Used different sized glassware joints to precisely control the flow of liquid
- Prevents raw liquid from moving from the flask to the condenser

Item Code	Cone Size to fit flask	Cone Size to fit Condenser	Pack of
030.202.01	14/23	14/23	1
030.202.02	19/26	24/29	1
030.202.03	29/32	29/32	1

ASTM Sizes

Item Code	Cap. (ml)	Socket Size	Cone Size	Pack of
030.402.04A	100	24/40	24/40	1

3.3
BORO
GLASS



1 x
SOCKET

1 x
CONE



ENQUIRE



ENQUIRE

