

Adapter, Receiver Delivery with short stem

- With 105° plain bend
- Used to extend drip tip to receiver flask

Item Code	Socket Size	Approx Stem Length (mm)	Pack of
019.202.01	14/23	65	10
019.202.02	19/26	65	10
019.202.03	24/29	65	10
019.202.04	29/32	65	10

3.3
BORO
GLASS



1 x
SOCKET



ENQUIRE



Adapter, Receiver Delivery with Long stem

- With 105° plain bend
- Used to extend drip tip to receiver flask

Item Code	Socket Size	Approx Stem Length (mm)	Pack of
019.202.05	14/23	190	10
019.202.06	19/26	200	10
019.202.07	24/29	200	10
019.202.08	29/32	200	10

3.3
BORO
GLASS



1 x
SOCKET



ENQUIRE



Adapter, Receiver, Straight

- Used to extend drip tip to receiver flask

Item Code	Socket Size	Pack of
020.202.01	14/23	10
020.202.02	19/26	10
020.202.03	24/29	10
020.202.04	29/32	10

3.3
BORO
GLASS



1 x
SOCKET



ENQUIRE

