

Flask, Boiling, Round Bottom, Single Neck with Joint

- Made from heat resistant ASTM E-438 Type 1, Boro 3.3 glass
- Complies with DIN / ISO 4797 & USP standards

Item Code	Capacity (ml)	Socket Size	Height (mm)	Approx Dia (mm)	Pack of
057.202.05	25	14/23	85	41	10
057.202.06*	25	19/26	85	41	10
057.202.07*	25	24/29	85	41	10
057.202.08	50	14/23	90	51	10
057.202.09	50	19/26	90	51	10
057.202.10	50	24/29	90	51	10
057.202.10A	50	29/32	90	51	10
057.202.11	100	14/23	105	64	10
057.202.12	100	19/26	105	64	10
057.202.13	100	24/29	105	64	10
057.202.14	100	29/32	105	64	10
057.202.16*	150	19/26	137	75	10
057.202.17*	150	24/29	137	75	10
057.202.18*	150	29/32	137	75	10
057.202.20*	250	14/23	140	85	10
057.202.21	250	19/26	140	85	10
057.202.22	250	24/29	140	85	10
057.202.23	250	29/32	140	85	10
057.202.24*	250	34/35	140	85	10
057.202.25*	500	19/26	163	105	10
057.202.26	500	24/29	163	105	10
057.202.27	500	29/32	163	105	10
057.202.28*	500	34/35	163	105	10
057.202.29*	1000	19/26	200	131	10
057.202.30	1000	24/29	200	131	10
057.202.31	1000	29/32	200	131	10
057.202.32*	1000	34/35	200	131	10
057.202.33	2000	24/29	240	165	6
057.202.34	2000	29/32	240	165	6
057.202.35*	2000	34/35	240	165	6
057.202.36A*	3000	24/29	260	185	1
057.202.36*	3000	29/32	260	185	1
057.202.37*	3000	45/40	260	185	1
057.202.40A*	5000	24/29	305	223	1
057.202.40*	5000	29/32	305	223	1
057.202.40B*	5000	34/35	305	223	1
057.202.41*	5000	45/40	305	223	1
057.202.44*	10000	29/32	380	279	1
057.202.45A*	10000	34/35	380	279	1
057.202.45	10000	45/40	380	279	1
057.202.45B*	10000	55/44	380	279	1
057.202.47B*	20000	55/44	435	345	1



3.3
BORO
GLASS

WHITE
PRINTING

Chemical
Resistant



ENQUIRE

*Non DIN