

# Condenser



## Condenser, Leibig with drip tip

- Complies with DIN 12576 standard
- Made From Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- With ground socket and cone joint
- With straight inner tube, outer joints and glass connector
- Resistant to chemicals and high temperatures

Item Code	Effective Length (mm)	Cone Size	Socket Size	Connector	Pack of
182.202.01	160	14/23	14/23	Glass	1
182.202.01A	160	19/26	19/26	Glass	1
182.204.02	200	19/26	19/26	Glass	1
182.202.02	250	19/26	19/26	Glass	1
182.202.03*	250	24/29	24/29	Glass	1
182.202.03A	250	29/32	29/32	Glass	1
182.202.04*	300	19/26	19/26	Glass	1
182.202.05*	300	24/29	24/29	Glass	1
182.202.06*	300	29/32	29/32	Glass	1
182.202.10*	400	24/29	24/29	Glass	1
182.202.11	400	29/32	29/32	Glass	1

## ASTM Size

Item Code	Effective Length (mm)	Cone Size	Socket Size	Connector	Pack of
182.402.01	120	14/20	14/20	Glass	1
182.402.02	180	19/22	19/22	Glass	1
182.402.04	200	24/40	24/40	Glass	1
182.402.05	300	24/40	24/40	Glass	1

3.3  
BORO  
GLASS

WHITE  
PRINTING

Chemical  
Resistant

1 x  
SOCKET

1 x  
CONE



ENQUIRE



ENQUIRE

\* Non DIN



## Condenser, Leibig with drip tip

- Complies with DIN 12576 standard
- Made From Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- With ground socket and cone joint
- With straight inner tube, outer joints and screw thread connector
- Resistant to chemicals and high temperatures

Item Code	Effective Length (mm)	Cone Size	Socket Size	Connector	Pack of
182.181.01	160	14/23	14/23	Screw Thread	1
182.181.01A	160	19/26	19/26	Screw Thread	1
182.181.02	250	19/26	19/26	Screw Thread	1
182.181.03*	250	24/29	24/29	Screw Thread	1
182.181.3A	250	29/32	29/32	Screw Thread	1
182.181.04*	300	19/26	19/26	Screw Thread	1
182.181.05*	300	24/29	24/29	Screw Thread	1
182.181.06*	300	29/32	29/32	Screw Thread	1
182.181.10*	400	24/29	24/29	Screw Thread	1
182.181.11	400	29/32	29/32	Screw Thread	1

## ASTM Size

Item Code	Effective Length (mm)	Cone Size	Socket Size	Connector	Pack of
182.481.01	120	14/20	14/20	Screw Thread	1
182.481.02	180	19/22	19/22	Screw Thread	1
182.481.03	200	24/40	24/40	Screw Thread	1
182.481.05	300	24/40	24/40	Screw Thread	1

3.3  
BORO  
GLASS

WHITE  
PRINTING

Chemical  
Resistant

SCREW  
THREAD



ENQUIRE



ENQUIRE

\* Non DIN